



Hearing Journey

Get Started

Need help? Call (779) 212 7334

This Quick Start Guide helps you set up your hearing aids. **Make sure to turn on your smartphone's Bluetooth®.**

Download the **FREE Lexie app**. Scan the QR code or find it on the App Store or Google Play.



STEP 1 Get to know your hearing aids




STEP 2 Power on and position

To power on, insert the battery, remove the sticker, and leave it open for 3-5 minutes before closing the battery door.

Place the hearing aid on your ear, wrap the wire around the front, and insert the dome.

Use the Lexie app to pair with your smartphone.



STEP 3 Selecting domes

Try different dome sizes to find what's most comfortable for you.



Open domes provide a natural sound.




Closed domes block external noise and may reduce feedback at high volumes.

STEP 4 Choosing receiver wires

Choose the size that fits comfortably and securely in your ear.

Use the included ear measuring tool to find the right size for each ear.



STEP 5 Adjust World Volume

World Volume changes how loud or quiet things sound. How to control it:

Via Lexie app: Drag the World Volume slider to the right to make things louder and to the left to make them quieter.

On hearing aid: Press the World Volume up or down buttons.

STEP 6 Customize your sound experience

Use the in-app adjustable features that allow you to adjust the World Volume of your two hearing aids, Treble/Bass, and Directionality. These options are available via the Controls screen of your Lexie app.

What To Expect

WEEK 1

Patience is key

It takes your brain and ears about 3-4 weeks to adjust to hearing aids.



Start with **1-2 hours** of wear time **once or twice daily** at low volume in quiet places. Then gradually increase volume and wear time.

Itchy ear canals or your own voice sounding different are normal and will likely fade over time.

WEEKS 2-3

Practice



Wear your hearing aids for **3-4 hours daily** in varied listening situations, from quiet to noisy, indoors, and outdoors.

Try different Environment settings based on your listening scenario.

WEEKS 3-4

Own your journey



Start increasing your wear time from **4 to 8 hours per day**.

Remember to take regular breaks while getting used to your hearing aids, especially if you're feeling overwhelmed.

WEEK 5 & BEYOND

Keep up the good work



Continue wearing your hearing aids **8+ hours daily** to reinforce your brain's learning.

Check your domes and wax guards every few days to see if they need cleaning.

Remove debris from the microphone inlets using the cleaning brush or a lint-free cloth. Consider using drying capsules to draw moisture from your devices, thereby prolonging their lifespan.

See "Maintain the hearing aids" in Section 8 of the User Manual for details.

Need help?

Call **(779) 212 7334**

Email **support@lexiehearing.com**

This document is only a brief overview to get you started. Please check the User Manual for more details. You can also reach us via phone or email for assistance.

Also note, if the pre-assembled size 2 receiver wire doesn't fit, please contact us directly to request the correct size. If you are within the 45-day trial period, you can receive a different size wire for free.

The Bluetooth® word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. and any use of such marks is under license. Other trademarks and trade names are those of their respective owners.

©2020-2025 hearX USA, Inc. All rights reserved. 44-0006-05, REV4

lexie®