

# Welcome to **lexie**

## Let's get started.

1



### Download the Lexie app.


Download from the App Store or Google Play Store. The app lets you control your hearing aids and teaches you how to wear them.

### Power on the hearing aids.

Press and release the bottom button on each of them.

**Note:** the app will confirm that they're powered on during setup.

2



### Follow the instructions in the app.

**Need help?** See reverse for more information.





**We're here to help.**

**Call (779) 212-7334**



For more detailed instructions, refer to the “Set up the hearing aids” section in the user manual.



Visit the Lexie product support webpage for videos, downloads, articles, and other resources.

**<https://lexiehearing.com/us/support>**

