

lexie

# Welcome

to Lexie Hearing



## Let's get started.



1

### Download the Lexie app.

Download from the App Store or Google Play Store. The app lets you control your hearing aids and teaches you how to wear them.



2

### Power on the hearing aids.

Press and release the bottom button on each of them.

Note: the app will confirm that they're powered on during setup.



3

### Follow the instructions in the app.

Need help? See reverse for more information.



## We're here to help

Call (800) 499-1336



For more detailed instructions, refer to section 2 in the user manual



Visit the Lexie product support webpage for videos, downloads, articles, and other resources.

<https://lexiehearing.com/us/support>

**POWERED BY**  
***BOSE***

©2020-2024 hearX IP (Pty) Ltd.

44-0011-03, REV2